


**BAUSCH+LOMB**  
**Instruments**  
Ear, Nose, Throat & Plastic Surgery

First Edition • US Only  
[www.BauschInstruments.com](http://www.BauschInstruments.com)

# Sinus Instruments Catalog



## Contents

Curettes .....	1
Elevators.....	1
Forceps	
Biopsy.....	2 - 3
Blakesley .....	3 - 5
Gruenwald .....	6
Heuwieser .....	7
Sinus.....	7 - 8
Takahashi.....	8
Knives .....	8
Needles .....	8
Punches	
Antrum Sidebiting .....	9
Circular Cutting .....	9
Ostrum Antrum .....	10
Rotating.....	11
Scissors .....	11
Specula.....	12
Suction Tubes .....	12

All original illustrations and text in this publication are fully protected by copyright, no part of it may be reproduced in any form without written permission by the publisher. Size specifications are stated as accurately as possible. Actual size may vary slightly due to handcrafted instrumentation.

©2012 Bausch & Lomb Incorporated. ®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates. 08/12.

# Sinus Instruments CURETTES, ELEVATORS

N2912

## Antrum Curette



- 1.5mm x 6mm cup.
- Overall length: 190mm, 7.5 inches.



N2750 1

## Coakley Antrum Curette Large Light Curve



- Large.
- Light Curve.
- 6mm x 8mm tip.
- Overall length: 174mm, 7.0 inches.



N2750 3

## Coakley Antrum Curette Large Strong Curve



- Large.
- Strong Curve.
- 6mm x 9mm tip.
- Overall length: 170mm, 6.8 inches.



N2750 5

## Coakley Antrum Curette Small Strong Curve



- Small.
- Strong Curve.
- 5mm x 7mm tip.
- Overall length: 168mm, 6.7 inches.



N2348

## Freer Septum Elevator 4.5mm



- Double-ended.
- Sharp and blunt 4.5mm tips.
- Length of shaft: 190mm, 7.5 inches.



### N2932 H **Biopsy Forceps 45° Horizontal 2 x 4mm**



- 45°, horizontal.
- Double action.
- Cup jaws, 2mm x 4mm.
- With Luer-Lock
- Working length: 130mm, 5.1 inches.

### N2932 V **Biopsy Forceps 45° Vertical 2 x 4mm**



- 45°, vertical.
- Double action.
- Cup jaws, 2mm x 4mm.
- With Luer-Lock
- Working length: 130mm, 5.1 inches.

### N2933 H **Biopsy Forceps 90° Horizontal 2 x 4mm**



- 90°, horizontal.
- Double action.
- Cup jaws, 2mm x 4mm.
- With Luer-Lock
- Working length: 130mm, 5.1 inches.

### N2933 V **Biopsy Forceps 90° Vertical 2 x 4mm**



- 90°, vertical.
- Double action.
- Cup jaws, 2mm x 4mm.
- With Luer-Lock
- Working length: 130mm, 5.1 inches.

### N2934 H **Biopsy Forceps 45° Horizontal 3 x 6mm**



- 45°, horizontal.
- Double action.
- Cup jaws, 3mm x 6mm.
- With Luer-Lock
- Working length: 130mm, 5.1 inches.

### N2934 V **Biopsy Forceps 45° Vertical 3 x 6mm**



- 45°, vertical.
- Double action.
- Cup jaws, 3mm x 6mm.
- With Luer-Lock
- Working length: 130mm, 5.1 inches.

### N2935 H **Biopsy Forceps 110° Horizontal 3 x 6mm**



- 110°, horizontal.
- Double action.
- Cup jaws, 3mm x 6mm.
- With Luer-Lock.
- Working length: 130mm, 5.1 inches.

### N2935 V **Biopsy Forceps 110° Vertical 3 x 6mm**



- 110°, vertical.
- Double action.
- Cup jaws, 3mm x 6mm.
- With Luer-Lock.
- Working length: 130mm, 5.1 inches.

### N2936 H **Biopsy Forceps 70° Horizontal 3 x 6mm**



- 70°, horizontal.
- Double action.
- Cup jaws, 3mm x 6mm.
- With Luer-Lock.
- Working length: 130mm, 5.1 inches.

### N2936 V **Biopsy Forceps 70° Vertical 3 x 6mm**



- 70°, vertical.
- Double action.
- Cup jaws, 3mm x 6mm.
- With Luer-Lock.
- Working length: 130mm, 5.1 inches.



# Sinus Instruments

## FORCEPS-BIOPSY, BLAKESLEY

N2931

### Biopsy Forceps Straight 2 x 4mm



- Straight.
- Double action.
- Cup jaws, 2mm x 4mm.
- With Luer-Lock
- Working length: 130mm, 5.1 inches.



N2990

### Blakesley Nasal Forceps No. 0



- No. 0.
- 3.5mm x 7mm cups, fenestrated.
- Length of shaft: 100mm, 3.9 inches.



N2991

### Blakesley Nasal Forceps No. 1



- No. 1.
- 4mm x 11mm cups, fenestrated.
- Length of shaft: 100mm, 3.9 inches.

N2993 1

### Blakesley Nasal Forceps No. 1 Upturned 45°



- No. 1.
- Upturned, 45°.
- 4mm x 10mm cups.
- Length of shaft: 105mm, 4.1 inches.

# Sinus Instruments

## FORCEPS-BLAKESLEY

### N3001 20 **Blakesley Through-Cutting Forceps Straight 2mm**



- Straight.
- 2mm bite.
- Length of shaft: 115mm, 4.5 inches.

### N3001 25 **Blakesley Through-Cutting Forceps Straight 2.5mm**



- Straight.
- 2.5mm bite.
- Length of shaft: 115mm, 4.5 inches.

### N3001 30 **Blakesley Through-Cutting Forceps Straight 3mm**



- Straight.
- 3mm bite.
- Length of shaft: 115mm, 4.5 inches.

### N3003 25 **Blakesley Through-Cutting Forceps Straight 2.5mm Serrated**



- Straight, serrated.
- 2.5mm bite.
- Length of shaft: 115mm, 4.5 inches.



### N3003 20 **Blakesley Through-Cutting Forceps Straight 2mm Serrated**



- Straight, serrated.
- 2mm bite.
- Length of shaft: 115mm, 4.5 inches.

### N3003 30 **Blakesley Through-Cutting Forceps Straight 3mm Serrated**



- Straight, serrated.
- 3mm bite.
- Length of shaft: 115mm, 4.5 inches.



# Sinus Instruments FORCEPS-BLAKESLEY

## N3002 20 **Blakesley Through-Cutting Forceps Upturned 2mm**



- Upturned.
- 2mm bite.
- Length of shaft: 115mm, 4.5 inches.

## N3002 25 **Blakesley Through-Cutting Forceps Upturned 2.5mm**



- Upturned.
- 2.5mm bite.
- Length of shaft: 115mm, 4.5 inches.

## N3002 30 **Blakesley Through-Cutting Forceps Upturned 3mm**



- Upturned.
- 3mm bite.
- Length of shaft: 115mm, 4.5 inches.

## N3004 25 **Blakesley Through-Cutting Forceps Upturned 2.5mm Serrated**



- Upturned, serrated.
- 2.5mm bite.
- Length of shaft: 115mm, 4.5 inches.



## N3004 20 **Blakesley Through-Cutting Forceps Upturned 2mm Serrated**



- Upturned, serrated.
- 2mm bite.
- Length of shaft: 115mm, 4.5 inches.

## N3004 30 **Blakesley Through-Cutting Forceps Upturned 3mm Serrated**



- Upturned, serrated.
- 3mm bite.
- Length of shaft: 115mm, 4.5 inches.

## N2897 **Blakesley-Wilde Ethmoid Forceps**



- 4mm x 10mm cups.
- Upturned.
- Length of shaft 110mm, 4.3 inches.



# Sinus Instruments FORCEPS-GRUENWALD

N9212

## Gruenwald Through-Cutting Forceps Curved 2mm



- Curved up.
- 2mm bite.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.



N9211

## Gruenwald Through-Cutting Forceps Curved 2.5mm



- Curved up.
- 2.5mm bite.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.

N9208

## Gruenwald Through-Cutting Forceps Curved 3mm



- Curved up.
- 3mm bite.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.

N9213

## Gruenwald Through-Cutting Forceps Straight 2mm



- Straight.
- 2mm bite.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.



N9210

## Gruenwald Through-Cutting Forceps Straight 2.5mm



- Straight.
- 2.5mm bite.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.

N9209

## Gruenwald Through-Cutting Forceps Straight 3mm



- Straight.
- 3mm bite.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.



N2905

### Heuwieser Antrum Grasping Forceps Curved



- Curved down.
- With Luer-Lock
- Working length: 130mm, 5.1 inches.



N2928 H

### Sinus Through-Cutting Forceps 45° Horizontal



- 45°, upturned.
- Horizontal.
- 2mm, double action jaws.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.



N2928 V

### Sinus Through-Cutting Forceps 45° Vertical



- 45°, upturned.
- Vertical.
- 2mm, double action jaws.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.

N2927 H

### Sinus Through-Cutting Forceps 90° Horizontal



- 90°, upturned.
- Horizontal.
- 2mm, double action jaws.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.

N2927 V

### Sinus Through-Cutting Forceps 90° Vertical



- 90°, upturned.
- Vertical.
- 2mm, double action jaws.
- With Luer-Lock.
- Length of shaft: 130mm, 5.1 inches.

N2898

### Sinus Cup Forceps



- 3mm cups.
- Length of shaft 180mm, 7.1 inches.



G6510

### Takahashi Nasal Forceps 2mm



- 2mm x 10mm jaws.
- Length of shaft: 162mm, 6.4 inches.



N2909

### Sickle Knife



- 9mm curved blade.
- Flat handle.
- Overall length: 190mm, 7.5 inches.



N2917

### Septum Needle



- Angled.
- 13 gauge shaft, 23 gauge tip.
- 14mm long needle tip.
- Overall length: 129mm, 5.0 inches.



# Sinus Instruments

## PUNCHES-ANTRUM, CIRCULAR

### N2903 L **Antrum Punch Sidebiting 180° Left**



- Through-cutting, left.
- 180°.
- 3.5mm bite
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.



### N2903 R **Antrum Punch Sidebiting 180° Right**



- Through-cutting, right.
- 180°.
- 3.5mm bite
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2904 L **Antrum Punch Sidebiting 40° Left**



- Through-cutting, left.
- 40°.
- 3.5mm bite
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2904 R **Antrum Punch Sidebiting 40° Right**




- Through-cutting, right.
- 40°
- 3.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2922 S **Circular Cutting Punch 3.5mm Straight**



- Straight, circular cutting
- 3.5mm bite.
- With Luer-Lock.
- Length of shaft: 160mm, 6.3 inches.



### N2922 U **Circular Cutting Punch 3.5mm Upturned**



- Upturned, circular cutting
- 3.5mm bite.
- With Luer-Lock.
- Length of shaft: 160mm, 6.3 inches.

### N2923 S **Circular Cutting Punch 4.5mm Straight**



- Straight, circular cutting
- 4.5mm bite.
- With Luer-Lock.
- Length of shaft: 160mm, 6.3 inches.

### N2923 U **Circular Cutting Punch 4.5mm Upturned**



- Upturned, circular cutting.
- 4.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2918 L **Ostrum Antrum Punch Pediatric Left 1.5mm**



- Backbiter, left side backward cutting.
- 1.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2918 R **Ostrum Antrum Punch Pediatric Right 1.5mm**



- Backbiter, right side backward cutting.
- 1.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2918 U **Ostrum Antrum Punch Pediatric Up 1.5mm**



- Backbiter, up side backward cutting.
- 1.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2919 L **Ostrum Antrum Punch Left 2.5mm**



- Backbiter, left side backward cutting.
- 2.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.



### N2919 U **Ostrum Antrum Punch Up 2.5mm**




- Backbiter, up side backward cutting.
- 2.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2919 R **Ostrum Antrum Punch Right 2.5mm**



- Backbiter, right side backward cutting.
- 2.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2896 1 **Ostrum Antrum Punch Forceps**



- Backward cutting.
- 3mm x 7mm bite.
- Length of shaft 100mm, 3.9 inches.



# Sinus Instruments

## PUNCHES-ROTATING, SCISSORS

### N2921 D Rotating Punch Small Curved Down



- Small, curved down, rotating.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2921 U Rotating Punch Small Curved Up



- Small, curved up, rotating.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2993 ROTP Rotating Punch Pediatric 1.5mm



- Pediatric, rotating.
- 1.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.



### N2994 ROT Rotating Punch Pediatric 2.5mm



- Pediatric, rotating.
- 2.5mm bite.
- With Luer-Lock.
- Length of shaft: 100mm, 3.9 inches.

### N2994 1 Nasal Sinus Scissors Straight



- Straight.
- 11mm blades.
- Length of shaft: 110mm, 4.3 inches.

### N2994 2 Nasal Sinus Scissors Left



- Left.
- 11mm blades.
- Length of shaft: 110mm, 4.3 inches.

### N2994 3 Nasal Sinus Scissors Right



- Right.
- 11mm blades.
- Length of shaft: 110mm, 4.3 inches.



N2101

### Nasal Speculum No. 1

- No. 1
- 30mm conventional blades
- Thin well-rounded tips.
- Overall length: 130mm, 5.1 inches.



N2906

### Suction Tube 3mm Long Curve

- Long curve.
- 3mm diameter.
- Bulbous tip.
- With thumb plate and cutoff (bypass hole).
- 25mm working length.
- Overall length: 135mm, 5.4 inches.



N2907

### Suction Tube 3mm Short Curve

- Short curve.
- 3mm diameter.
- Bulbous tip.
- With thumb plate and cutoff (bypass hole).
- 13mm working length.
- Overall length: 146mm, 5.7 inches.

N2908

### Suction Tube 4mm Short Curve

- Short curve.
- 4mm diameter.
- Bulbous tip.
- With thumb plate and cutoff (bypass hole).
- 10mm working length.
- Overall length: 146mm, 5.7 inches.



Antrum Curette.....	1	N2922 S.....	1, 9
Antrum Punch Sidebiting 40° Left.....	9	N2922 U.....	9
Antrum Punch Sidebiting 40° Right.....	9	N2923 S.....	9
Antrum Punch Sidebiting 180° Left.....	9	N2923 U.....	9
Antrum Punch Sidebiting 180° Right.....	9	N2927 H.....	7
Biopsy Forceps 45° Horizontal 2 x 4.....	2	N2927 V.....	7
Biopsy Forceps 45° Horizontal 3 x 6.....	2	N2928 H.....	7
Biopsy Forceps 45° Vertical 2 x 4.....	2	N2928 V.....	7
Biopsy Forceps 45° Vertical 3 x 6.....	2	N2931.....	3, 5, 8
Biopsy Forceps 70° Horizontal 3 x 6.....	2	N2932 H.....	2
Biopsy Forceps 70° Vertical 3 x 6.....	2	N2932 V.....	2
Biopsy Forceps 90° Horizontal 2 x 4.....	2	N2933 H.....	2
Biopsy Forceps 90° Vertical 2 x 4.....	2	N2933 V.....	2
Biopsy Forceps 110° Horizontal 3 x 6.....	2	N2934 H.....	2
Biopsy Forceps 110° Vertical 3 x 6.....	2	N2934 V.....	2
Biopsy Forceps Straight 2 x 4.....	3, 5, 8	N2935 H.....	2
Blakesley Nasal Forceps No. 0.....	3	N2936 H.....	2
Blakesley Nasal Forceps No. 1.....	3	N2936 V.....	2
Blakesley Nasal Forceps No. 1 Upturned 45°.....	3	N2990.....	3
Blakesley Through-Cutting Forceps Straight 2.5mm.....	4	N2991.....	3
Blakesley Through-Cutting Forceps Straight 2mm.....	3, 4	N2993 1.....	3
Blakesley Through-Cutting Forceps Straight 3mm.....	4	N2993 ROTP.....	11
Blakesley Through-Cutting Forceps Upturned 2.5mm.....	5	N2994 1.....	11
Blakesley Through-Cutting Forceps Upturned 2mm.....	5	N2994 2.....	11
Blakesley Through-Cutting Forceps Upturned 3mm.....	5	N2994 3.....	11
Blakesley-Wilde Ethmoid Forceps.....	5	N2994 ROT.....	11
Circular Cutting Punch 3.5mm Straight.....	1, 9	N3001 20.....	3, 4
Circular Cutting Punch 3.5mm Upturned.....	9	N3001 25.....	4
Circular Cutting Punch 4.5mm Straight.....	9	N3001 30.....	4
Circular Cutting Punch 4.5mm Upturned.....	9	N3002 20.....	5
Coakley Antrum CuretteLarge Light Curve.....	1	N3002 25.....	5
Coakley Antrum CuretteLarge Strong Curve.....	1	N3002 30.....	5
Coakley Antrum Curette Small Strong Curve.....	1	N3003 20.....	4
Freer Septum Elevator 4.5mm.....	1	N3003 25.....	4
G6510.....	8	N3003 30.....	4
Gruenwald Through-Cutting Forceps Curved 2.5mm.....	6	N3004 20.....	5
Gruenwald Through-Cutting Forceps Curved 2mm.....	6	N9208.....	6
Gruenwald Through-Cutting Forceps Curved 3mm.....	6	N9209.....	6
Gruenwald Through-Cutting Forceps Straight 2.5mm.....	6	N9210.....	6
Gruenwald Through-Cutting Forceps Straight 2mm.....	6	N9211.....	6
Gruenwald Through-Cutting Forceps Straight 3mm.....	6	N9212.....	6
Heuwieser Antrum Grasping Forceps Curved.....	7	N9213.....	6
N2101.....	12	Nasal Sinus Scissors Left.....	11
N2348.....	1	Nasal Sinus Scissors Right.....	11
N2750 1.....	1	Nasal Sinus Scissors Straight.....	11
N2750 3.....	1	Nasal Speculum No. 1.....	12
N2750 5.....	1	Ostrum Antrum Punch Forceps.....	10
N2896 1.....	10	Ostrum Antrum Punch Left 2.5mm.....	10
N2897.....	5	Ostrum Antrum Punch Pediatric Left 1.5mm.....	10
N2898.....	8	Ostrum Antrum Punch Pediatric Right 1.5mm.....	10
N2903 L.....	9	Ostrum Antrum Punch Pediatric Up 1.5mm.....	10
N2903 R.....	9	Ostrum Antrum Punch Right 2.5mm.....	10
N2904 L.....	9	Ostrum Antrum Punch Up 2.5mm.....	10
N2904 R.....	9	Rotating Punch Pediatric 1.5mm.....	11
N2905.....	7	Rotating Punch Pediatric 2.5mm.....	11
N2906.....	12	Rotating Punch Small Curved Down.....	1, 11, 12
N2907.....	12	Rotating Punch Small Curved Up.....	11
N2908.....	12	Septum Needle.....	8
N2909.....	8	Sickle Knife.....	8
N2912.....	1	Sinus Cup Forceps.....	8
N2917.....	8	Sinus Through-Cutting Forceps 45° Horizontal.....	7
N2918 L.....	10	Sinus Through-Cutting Forceps 45° Vertical.....	7
N2918 R.....	10	Sinus Through-Cutting Forceps 90° Horizontal.....	7
N2918 U.....	10	Sinus Through-Cutting Forceps 90° Vertical.....	7
N2919 L.....	10	Suction Tube 3mm Long Curve.....	12
N2919 R.....	10	Suction Tube 3mm Short Curve.....	12
N2919 U.....	10	Suction Tube 4mm Short Curve.....	12
N2921 D.....	1, 11, 12	Takahashi Nasal Forceps 2mm.....	8
N2921 U.....	11		

# **BAUSCH + LOMB**

# **Instruments**

Ear, Nose, Throat & Plastic Surgery



To request a catalog of our full line of Ear, Nose, Throat & Plastic Surgery Instruments, visit [www.bauschinstruments.com](http://www.bauschinstruments.com) or contact Customer Service at (800) 338-2020.

For the most up-to-date listing of available instruments, visit [www.bauschinstruments.com](http://www.bauschinstruments.com)

**Customer Service:**  
**(800) 338-2020**

**Fax: (800) 362-7006**

**International Fax: (636) 226-3119**

**Instrument Repairs:**  
**(800) 658-4751**

**Custom Instruments:**  
**(800) 325-9505 x3779**

All original illustrations and text in this publication are fully protected by copyright, no part of it may be reproduced in any form without written permission by the publisher. Size specifications are stated as accurately as possible. Actual size may vary slightly due to handcrafted instrumentation.

©2012 Bausch & Lomb Incorporated. ®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates. 08/12. SU6795